NEWYES[®] & EYECOS[®]

Eye Aesthetics

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Estética Ocular

Pioneers of the Future

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Pioneros del Futuro

New laser technology to change eye color
Eye Aesthetics

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Summary
Introduction
Material and methods
Classification of pigmentation levels
Results
Photographic Series
Complications
Discussion
Conclusions

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SUMMARY

RESUMEN

Neweyes is the first laser technique capable of changing the color of the iris and, therefore, of the eyes, in a safe, effective and predictable manner. The results obtained are all natural, and generally predictable.

Neweyes is a non-surgical outpatient laser technique, and thus free of typical intraocular surgery complications.

The most commonly observed practical complication is the stromal micro bleeding, although it is self-limited to 100% of the cases. It resolves spontaneously, while leaving no alterations of any kind.

Ocular hypertension is the most serious complication in the short, medium and long-term. The present review outlines the measures to prevent it, minimize and manage it almost completely, avoiding the appearance of anatomical and functional sequels.

In a series of 180 eyes treated, patients have reported no significant complaints of discomfort, pain, irritation, dryness, photophobia or decreased vision. There have been no elevations in intraocular pressure neither refractive changes.

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INTRODUCTION

Humans greatly appreciate originality, uniqueness, singularity, rarity and exoticism, as goods prized for their scarcity and the meaning of power that they involve.

The eyes, the first traits explored when meeting a person, are a powerful tool for communication, for seduction and even as a power source. And the rarer the color, the more effective their subliminal message. The least common eyes are the result of a specific genetic trend or even of a disease.

The eyes are called "mirrors of the soul", and although we do not realize it, they are continually used to express emotions, moods and even to make others to fall in love.

The eye color is a genetic trait that is determined by the amount and distribution of melanin in the iris. There are three elements that contribute to give iris coloration: posterior epithelium melanin, the stroma melanin and anterior epithelium melanin.

Besides melanin, lipochrome also acts in the process of pigmentation. There are two types of melanin, eumelanin, with dark brown appearance, and pheomelanin, which is reddish-yellowish.

Posterior epithelium pigment is always eumelanin, present in large quantities. Therefore, the color of the eyes is determined by the anterior epithelium pigment and the stromal tissue density.

In the world there is a relative uniformity in eye color, brown being the most frequent. However, the exception is Europe, where a diversity have been generated that encompasses the hazel, green, blue and gray.

Reasons were made of natural selection, hybridization with Neanderthals men, and even sexual selection (Luigi Cavalli-Sforza), according to which the brightest and gaudiest colors would have involved advantages when choosing a mate.

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GENETIC DETERMINANTS

The eye color is an inherited character influenced by several genes:

- EYCL1: On chromosome 19, responsible for green and blue, dominant and recessive allele.
- EYCL2: On chromosome 15, responsible for the brown.
- EYCL3: Regulates the amount of melanin in the body.
- OCA1-4: They are located on chromosome 15. They are responsible for the pigmentation of skin and hair, and when mutated or altered they cause various types of ocular skin albinisms. The OCA2 induces brown and blue tones, dominant and recessive allele.
- SLC24A4.
- TYR.

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TYPES OF COLORS IN HUMANS

→ The tan, brown or chestnut contains large amounts of melanin in the iris frontal epithelium, to absorb a large amount of short-wavelength light.

→ Amber is a mixture of brown and hazel. It is generated by the lipochrome.


→ Hazelnut color is intermediate between brown and green.

→ Green is mixture of brown and blue, or grey. In green eyes there are moderate amounts of melanin.


→ Blue has very low densities of melanin pigment in the anterior epithelium. The stroma, which is a transparent tissue composed of collagen fibers, reflects outwards the light incident on the posterior epithelium. The light rays split across the stromal fibers and reflect the blue. Blue eyes first appeared 6,000 to 10,000 years in the northeast of the Black Sea, as a result of the mutation of a single individual. At present they are relatively common in Europe and other areas of continents with predominantly European ancestry.

The world population estimated at about 7,000 billion, has some 150 million people with blue eyes. That is, about 2% of the whole.

→ Grey. Caused by a high density of stromal tissue.

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BIRTH

Newborns often have depigmented eyes with a bluish appearance, but as they receive sunlight, the iris melanocytes begin to produce melanin, and the color switches to the adult coloring tone, reaching 50% of its intensity at 3 years old.

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→ Heterochromia : It is the case when both irises have different colors (full), or when in a single iris there are areas of different colors (partial or sectorial). It can be congenital or acquired. Congenital may be associated with Neurofibromatosis, Waardenburg syndrome or Claude Bernard - Horner syndrome. It is common in dogs, cats and horses.

→ Aniridia: bilateral, heritable and unusual. It is the total absence of the iris caused by an alteration of chromosome 11. The affected persons have low vision and photophobia, cataract, glaucoma, nystagmus, and macular and foveal hypoplasia.

→ Albinism: It is the complete absence of melanin in the skin, eyes and hair. In the severe grade, the eyes have red or even purple appearance. Affected persons suffer from photophobia, nystagmus, strabismus, myopia and astigmatism.

→ Kayser- Fleischer ring: Consists of copper deposition in the corneal periphery, and it has golden-greenish color. It is part of the Wilson 's disease or hepato - lenticular degeneration (widespread deposition of copper in tissues: liver, nervous system, cornea,...) and it is present with psychiatric symptoms, given the high toxicity to the tissue of the mineral.

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MEDICAL-OCULAR IMPLICATIONS

- ➔ A study by the University of Otago in New Zealand, has established the relationship between pigmented eyes and a greater predisposition to glaucoma.
- ➔ A clinical review of Eye & Ear Infirmary of Massachusetts, has found higher rates of deaths caused by systemic metastases in individuals with low pigmented eyes with choroidal melanoma.
- ➔ In the University of Detroit a higher prevalence of cases of AMD (age related macular degeneration) in cases of low iris pigmentation.

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METHODS FOR CHANGING THE EYE COLOR

→ Photoshop Software: One of the simplest yet most effective ways to transform the image of a person is changing the color of its eyes using Photoshop, in only 3 minutes, but only at photography level.

→ Cosmetic contact lenses: The cosmetic or decorative contact lenses change the color appearance of the eyes, and in case they are circular they make the iris to appear larger too.

→ Biokinesis: It is defined as the manipulation of the body through the mind. It consists of the use, repeated and prolonged in time, of auto-suggestion techniques with virtual images, which eventually causes to believe that there has been an actual physical modification..

→ Newiris: The process "New Color Iris" was developed by Dr. Kahn in Panama in 2002 to help treating the Kuna Indians, who have a high prevalence of albinism. During surgery, a 9-10 mm diameter implant is inserted in the eye, to coat the original iris.

→ Jewelry Implant in the Eye: The Institute for Innovative Ocular Surgery in Rotterdam performed the first 10 implants and published them in 2003 in the Journal of Cataract & Refractive Surgery which showed that the placement of platinum jewellery in the ocular globe generated no complications months later. The method, called Safe Sight Eye Jewelry, consists of performing, after local anesthesia of the eye, a small incision in the conjunctiva. Subsequently, a thin piece of jewellery is inserted into the incision. In total, it takes about 10-15 minutes.

→ Stroma Medical: The California based company Stroma Medical has been announcing from several years a laser able to change the eye color, irreversibly. The new treatment could be available in just three years after being approved by the FDA (the U.S. medical authority). To this end, should also raise the \$ 750,000 he needs to launch the production of this technique, which is called "Lumineyes" and which would be able to transform brown eyes to blue eyes. To this end, the application of several sessions lasting 20 seconds over a month would suffice.

→ Bright Ocular: A pending U.S. patent application, relating to an colored intraocular implant, thin, flexible, biocompatible, developed to alter iris appearance for medical and cosmetic purposes. It treats ocular albinism and heterochromia protecting the eyes from the harmful effects of sunlight. This artificial iris implant completely hiddens the iris defects such as coloboma, severe atrophy, iridoschisis, and is useful for patients whose only desire is

to change the color of their eyes for aesthetic reasons. The implants are made of silicone manufactured by Nusil, global biotechnology leader based in Santa Barbara, California. The method is outpatient, is performed under topical anesthesia and takes about 15 minutes for each eye. The implant can be inserted through a corneal incision of 2.8 mm, without sutures.

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MATERIAL AND METHODS

The NEWEYES technique has been applied to 90 patients, 49 women (54.4 %) and 41 men (45.5 %). A total of 180 eyes were treated.

The overall mean age was = 33.47 years. The average age of women was =33.8 years, and men = 32.65 years.

Regarding refractive errors, 54 patients (60%) were emmetropic (+ / - 0.50 diopters), myopic were 20 (22.2%), 12 had astigmatism (13.3%), and only 4 were long-sighted (4.4%).

The average eye pressure before treatment was = 12mmHg. Women exhibited = 11.9mmHg average pressure, and men = 12.1mmHg.

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PROTOCOL FOR CANDIDATES ADMISSION

- All those interested in receiving the NEWEYES laser that suffered any chronic general disease requiring continuous treatment, such as diabetes, collagen diseases, heart diseases, kidney diseases, immune deficiencies, depression, psychopathy,... have been discarded
- All those interested in receiving the NEWEYES laser that suffered any chronic eye disease requiring continuous treatment, such as glaucoma, uveitis, cataracts, keratopathies, retinopathy, neuropathy,... have been discarded.
- All those interested persons with direct family history with any chronic eye disease requiring continuous treatment, such as glaucoma, uveitis, cataracts, keratopathies, retinopathy, neuropathy,... have been discarded.
- Patients candidates for NEWEYES with refractive errors: myopia, long-sightedness, astigmatism and presbyopia, including two cases of partial amblyopia with anisometropia have been accepted.
- Only patients with an eye pressure lower than 17mmHg, measured with Goldmann tonometer, have been accepted as candidates for the study.
- Finally, only patients belonging to levels I, II and III of the Classification of levels of pigmentation (see section 4) have been accepted. All those interested with a level IV (excessive amount of pigment or too dark coloring) have been discarded.

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PRELIMINARY PROTOCOL FOR NEWEYES TREATMENT

1 The first day the patient is informed about the nature of the NEWEYES laser treatment.

2 It is provided with the Informed consent form that has to be read in the following days such that he can ponder the pros and cons, as well as consulting with family and friends. This document and other specific information is public and can be downloaded from the EYECOS website: "www.neweyeslaser.com"

3 High resolution photographs of both eyes are taken and shown to the patient concerned to explain the particulars of its iris, as regards color, pigment density and surface irregularity.

4 A prediction of the result is performed after completion of the treatment, based on the comparative extrapolation with cases of patients already treated (previous, on course and final photographs). It is advised that color, hue or regularity CANNOT BE CHOSEN.

5 Finally, the patient is clearly advised about the irreversibility of the treatment and that there are little experience data on it (24 months = 2 years), thus being impossible to be certain of its long term total safety.

6 The second day a previous basic eye examination is performed, to confirm that the candidate fulfills all the established conditions of admission.

7 The signed Informed Consent is obtained from the patient, thus obtaining the patient permission to start the laser sessions.

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LASER SESSIONS PROTOCOL

- 1** The first day 1 tablet of ibuprofen is administered and 2% pilocarpine drops are instilled in both eyes one hour before application.
- 2** Just before, prescaina anesthetic drops are instilled in both eyes, and the patient is instructed to cooperate to maintain the maximum fixation and immobility.
- 3** A “peripheral gonioplasty” is performed in both eyes to maximize the amplitude of the camerular angle.
- 4** No post laser treatment is administered.
- 5** The second day of laser and one hour before application, 1 tablet of ibuprofen is administered and 2% pilocarpine drops are instilled in the right eye.
- 6** Just before, prescaina anesthetic drops are instilled and an “iridoplasty” is carried out on $\frac{3}{4}$ middle portions of the surface of the iris, leaving untreated the peri-pupillary central and peripheral areas.
- 7** As post laser treatment two drops of Iopimax (0.50 %apraclonidine) are instilled, and scheduled as follows: Combigan eye drops (brimonidine + timolol: 1 drop / every 12 hours) and Nevanac eye drops (nepafenac: 1 drop / every 8 hours) for 1 week.
- 8** It is recommended to take for 48 hours, 1 ibuprofen tablet every 8-12 hours.
- 9** The patient is advised to call urgently if it has intense eye pain that does not relief after the first hours, or if progressive dilation of the pupil occurs. Both signs indicate a sharp rise in the intraocular pressure.
- 10** The third day of laser, and prior to the implementation thereof, the patient is asked about possible discomfort and his visual acuity, his biomicroscopy and his ocular pressure are explored. Photographs of the partial results, which are displayed to the person concerned will also be taken.
- 11** Therefore an “iridoplasty” is carried out in $\frac{3}{4}$ medial parts of the iris of the left eye. The pre and post treatment are the same as in the right eye.
- 12** The third and fourth laser days “iridoplasties” of the peripheral and remaining peri-pupillary areas are carried out, first the right and then the left eye. The pre and post treatments are identical.
- 13** After completing the first 5 NEWEYES laser sessions, a treatment is prescribed consisting of Hylo Parin eye drops (sodium hyaluronate + sodium heparin: 1 drop / every 8 hours), and the patient is appointed for review after 1 month, provided there is no discomfort. The patient may have a normal life, as regards physical and occupational or sports. It is not necessary to use protective sunglasses.
- 14** One month after the initial sessions, a visit with ocular exploration and shooting is done. If necessary, a final session of “point iridoplasty” is applied, to remove traces of residual pigment. The common pre and post treatment are applied.
- 15** From this point, the patient is treated several months with Hylo Parin eye drops, every 8 hours a day, and routine monitoring visits are conducted every month.

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CLASSIFICATION OF PIGMENTATION LEVELS

CIn order to properly classify the type of pigmentation of the iris of patients interested in receiving the NEWEYES laser, and to assess their level of intrinsic difficulty, a classification of pigmentation levels has been designed.

For setting the levels, three factors were taken into account:

- 1** Color Intensity: light, medium and dark brown.
- 2** Thickness of pigment: superfine, fine, medium and coarse.
- 3** Surface irregularity Level: none, medium and extreme.

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THE PIGMENTATION LEVELS ARE THE FOLLOWING

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I Bright Iris, with superfine or thin layer, and no irregularity.

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


II Bright or intermediate iris, with superfine or thin layer, and no irregularity.

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
III Intermediate iris, with medium or thick layer, and medium irregularity.

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IV Very dark iris or extreme irregularity.

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RESULTS

The parameters that define the quality and guarantee level of a new laser or surgical technique are:

1 Security **2** Effectiveness **3** Predictability

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

SECURITY

The NEWEYES laser of EYECOS is highly safe, as no refractive changes, decreased visual acuity, limitations in pupillary reflexes, subacute or chronic uveitis, or intraocular pressure elevations have been observed after the treatment has been completed.

While there is a blurring of vision for approximately 4-5 hours following each laser session, owing to the miotic effect of pilocarpine and the dispersion of melanin pigment in the anterior chamber, the vision is fully recovered without remarkable inconveniences. On three occasions, for persisting a soft iritis with traces of pigment in the aqueous humor, it took up 48 hours to recover vision.

Regarding the nerve – muscle functionality of the iris treated with NEWEYES, it has to be pointed out that it is fully recovered quickly after each session. Only in two patients a very soft unilateral midriasis showed up, with discreet anisocoria, which normalized after one week. Pupillary reflexes to light, both consensuated direct and indirect, remained unchanged after application of the laser.

Sunlight has not been a problem for patients. In general no complaints were made about it except on two occasions where a strongly increased sensitivity appeared owing to intense sunny days. In no case a frequent need to wear sunglasses after completion of treatment has been reported.

As for physical and occupational activity, its restriction is recommended only for 4-5 hours immediately following the laser sessions, going back to fully normal life the day after.

With respect to intraocular pressure, moderate tension peaks were observed in the first 3-4 hours after the completion of the first two sessions of “iridoplasties” because in these sessions larger iris areas are treated (3 / 4), and therefore, they are accompanied with an increased release of pigment. No elevation peaks were recorded after the completion of the “gonioplasties” or peripheral “iridoplasties”, peripupillary or finishing touches.

While the overall average previous pressure was = 12.1mmHg, the tension peaks had elevations about 10 to 15 mmHg, with a mean = 25.6mmHg. These elevations are similar to those found after performing trabeculoplasty with argon laser, so we have established as a protocol a similar preventive treatment based on Iopimax eye drops (0.50 %apraclonidine), immediately after the session, and Combigan and Nevanac eye drops, over one week. Ibuprofen tablets (every 8 - 12h) also have been added, for 48 hours so as to reduce the iris and frontal chamber inflammation.

After finishing a week of treatment with Combigan and Nevanac eye drops, these were completely stopped, the patient just having to be treated with Hylo Parin eyedrops (1 drop / every 8 hours) for 3-6 months.

In the short, medium and long term, and without topical hypotensive treatment, no intraocular pressure elevations taken with the Goldman tonometer have been recorded. The overall average, once the session completed, has been = 12.0mmHg, without significantly statistical differences with the pre-treatment pressures(12.1mmHg).

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NEWYES EFFECTIVENESS

The level of effectiveness of NEWYES laser is very high. Instantly a change in iris color is achieved, from brown to green or blue.

In 100% of the cases the objective sought has been achieved in completely removing the melanin pigment in the anterior iris, so that the eye color of the patient has been obviously apparently transformed, although discrete, from the initial color. In all cases the obtained hue was visibly clearer than the original brown.

On average, two NEWYES laser sessions per eye sufficed to achieve greenish or bluish colors. Usually, a bilateral review session has been performed after one month, so as to finish the removal of residual pigment traces, which are usually located in the peri pupillary area.

The power output from the NEWYES laser was proportional to the level of pigmentation of the iris. Very low in eyes with pigmentation level I and II, and moderate in level III. The minimum energy capable of achieving the iris depigmentation has been always used.

In the same eye, due to the irregularities in the distribution of the pigment, different power levels were used, generally higher on the area surrounding the pupil and lower in peripheral sectors.

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EFFECT AFTER THE FIRST NEWYES SESSION

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EFFECT AFTER THE SECOND NEWYES SESSION

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
EFFECT AFTER COMPLETING THE NEWYES TREATMENT

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NEWYES PREDICTABILITY

Having already demonstrated the safety and efficacy of the NEWYES laser, we must address the issue of the predictability level.

Before starting the treatment the patient is informed that he CAN NOT CHOOSE the final color. This color can vary from green to blue, with their respective range of hues.

While the practical NEWYES predictability is 100 % in achieving the pigmentation of the iris, is not as high with regards to the prediction of the exact color and tone at the end of treatment.

It is us, not the patient, who will inform, as accurately as possible, the likely outcome to be achieved. To do this we rely on three prediction factors:

1 Pigmentation level of the patient.

The lower the level, I and II, the more likely we are to get very bright green and blue colors. The higher the pigmentation, the more frequent green colors with dark shades will be.

2 Photographs comparison with other cases already made.

Photographs will be shown to the patient similar to him, with the same pigmentation level, so he will be able to check out the final result and have an idea of its possibilities.

3 Genetic characteristics of relatives.

As a general rule, where the family have green eyes the patient will end with green eyes. Conversely, if blue eyes abound among immediate family there is a high chance of achieving blue eyes.

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EXAMPLE Of a brown eye that turns to blue after neweyes. It was a case of congenital heterochromia, wherein the right was blue and the eye left brown (III)

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
EXAMPLE Of level I eye with uniform final result

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



EXAMPLE Examples of level II cases with bright green eyes results

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

EXAMPLE Case examples with light blue eyes results

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

EXAMPLE Case examples with bright blue eyes results

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

EXAMPLE Case examples with greenish dark grayish eyes results

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EXAMPLE Examples of looks with green and light blue eyes results

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








COMPLICATIONS

The NEWEYES security level is very high. firstly, there appears none of the typical complications of surgical intervention, such as infections, wound healing abnormalities, or relating to the suture.

Secondly, no adverse reactions to implants, prosthesis or chemicals and drugs occur, since none of them is introduced into the eyeball.

Only a laser is applied on the surface of the iris, with low power radiation without damaging any intraocular structure serving to vision, such as the cornea, lens, vitreous, retina, or the optic nerve.

Among the 180 eyes treated, we have observed two types of complications:

- 1** Those related to ocular hypertension.
- 2** Those derived from the direct action of the laser on the iris.

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OCULAR HYPERTENSION

The elevation of the acute and chronic pressure must be controlled. There is a tensional peak the first 3-4 hours after each session, which can be estimated on average between 10 and 15 mmHg. Our protocol is similar to that used in trabeculoplasty with argon laser, that is, with early instillation with Iopimax (0.50 % apraclonidine), use of Combigan (topical hypotensors) and Nevanac (topical anti-inflammatory), a week, and ingestion of ibuprofen tablets for two days. The patient must be advised to go to casualty department in case of long lasting severe headache, accompanied by non-reactive mydriasis, both signs of post laser acute glaucoma. In this case, the appropriate treatment has to be established without delay for quick resolution to avoid visual loss and associated damages.

Subacute and chronic hypertension has never been detected in routine checks, although the hypotensive patients stop the treatment just after one week. Undoubtedly, the long-term monitoring of pressure is necessary to foresee its possible elevation.

During the last fifteen years the anti glaucomatous therapeutic has greatly advanced, and thus the rates of filtering surgery have decreased, given the extraordinary effectiveness of the new generation of drugs, such as Latanoprost.

Surely we have enough tools to control pressure that may appear with age in patients treated with NEWEYES laser: topical medication, argon laser trabeculoplasty and modern depot implants. The deep non-perforating esclerectomy has also reduced the incidence of complications of classic perforating trabeculectomy, already in disuse.

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COMPLICATIONS ARISING FROM THE DIRECT ACTION OF THE LASER ON THE IRIS

This type of complication is relatively common, although it lacks severity and clinical relevance. It consists of setbacks that auto resolve spontaneously, leaving no functional sequels. They are the following:

- Stromal micro bleeding. They auto resolve within a few seconds with the help of ocular pressure.
- Lower pigment deposition (“6h goop”). Disappears over time, or a minimum rest remains attached to the lower iris, just 6 hours after.
- Perforating microiridotomy. Without transcendence.
- Mild anisocoria, caused by a low unilateral pupillary dilation. It normalizes in a few days or weeks.
- Blurred vision for 4-5 hours.
- Temporary photophobia for 1 week..
- Anterior uveitis, which resolves in 24, 48 or 72 hours.
- Urrets - Zabalía syndrome.

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

1 JUSTIFICATION

The huge and inexhaustible potential market that exists worldwide for lightening the color of the eyes, that is, changing the iris brown to greenish or blueish hues, does not justify the use of any technique for this purpose.

The urgency to achieve the magical transformation that interested patients “are begging for” should not be an excuse to expose them to surgeries with dubious effectiveness and that are highly dangerous.

Candidates argue family, professional, psychological, personal reasons... and must be the responsibility of the physician to differentiate between the legitimate and justified, realistic and grounded, safely achievable aspirations and the cases lacking of logic, in which fantasy and fictional imagination predominates over prudent rationality.

Thus, a judicious selection is essential. We have to discriminate accurately the interested patients and not accept anyone who wants to perform the technique in an impulsive manner and without enough thought.

The patients with systemic diseases who are subjected to chronic treatments will be discarded, because these treatments may interfere with the clinical course of the NEWEYES sessions thus favoring intraocular inflammation and the emergence of other complications. Those patients that exhibit prior medical pathologies, both personal and immediate family, type of glaucoma, retinopathy, uveitis, neuropathy,... will be discarded as well. On the contrary, the different types of graduation defects such as nearsightedness, farsightedness, astigmatism and presbyopia and those already undergoing refractive surgery such as Lasik will be accepted.

While the desire to change the personal image may seem legitimate, our duty is to advise and guide properly following the dictates of good medical practice, about the real results targeted, such as not to deceive the aspirations of patients, so many years treasured in the depths of their illusions.

At ophthalmologic level, we should not accept eyes with excessive pigment, because the treatment would increase the likelihood of complications such as uveitis and pigmentary glaucoma. And of course, we should also discourage the technique in cases of absence of iris pigment, in which the person intends to brighten his eyes irrationally.

Psychological disorders are a serious problem for the NEWEYES treatment. Anxiety, depression and obsessive perfectionism are not good allies. We contraindicate the technique in the presence of schizoid psychopathy, deep endogenous depression, expansive mania, and bipolar syndromes. These patients distort reality, and will never be satisfied with the results. Additionally, they could suffer a depressive crisis as a result of receiving the treatment, due to an acute decompensation of their mental illness.

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2 VERY DANGEROUS AND UNRELIABLE TECHNIQUES

Firstly, let's mention cosmetic anterior chamber Intra Ocular Lenses. Two are being commercialised currently on the market, although their use is very restrictive.

Dr. Kahn developed in 2002 in Panama, where he continues to make interventions, the New color iris project. It involves implanting foldable silicone lens with uniform coloring, with three main hues.

While it may seem a simple intervention, which is done through a small incision of 2.8-3 mm and with topical anesthesia, complications begin soon in the immediate postoperative with inflammation, visual discomfort and pain.

By introducing a large lens in the anterior chamber, a constant rubbing against the cornea and iris is generated, and a serious conflict of space is produced, which has an impossible resolution, unless the lens is extracted. In fact, there have been numerous reported cases of emergency implant extractions, early or late.

It is easy to imagine that the Newiris lenses will cause acute and chronic elevation of intraocular pressure by blocking the camerular angle and inflammatory uveitis with frontal synechiae. We do not know the frequency of secondary cataract, but it seems to us that it may be high.

From the aesthetic and functional point of view, the Newiris lenses have two distinct disadvantages. Firstly, the artificial and "jazzy" color

of the lenses, which in no way imitates the natural surface of the iris. Secondly, the entire limitation miosis - mydriasis gap of a normal pupil, in response to different lighting conditions. The lenses have a fixed pupillary aperture of 3-4 mm, so they never dilate in darkness or constrict in presence of sunlight. Therefore, the patient will lack sufficient light at night, and will be dazzled during bright days. Moreover, the visual field is permanently limited, and therefore, the driving and walking will suffer greatly.

The company Bright Ocular, headquartered in California, produces very similar lenses to Newiris, constructed of silicone for ophthalmic use, but they just offer a simple improvement, a slightly more natural appearance (similar to contact lenses) with a greater range of colors. Therefore, they present the same risk of serious complications than Newiris implants and we completely discourage their use.

Regarding Stroma Medical company, based also in California, and who for years much publicity has made of its Lumineyes laser, we are still waiting for a single photograph of his "amazing and Hollywood" results.

According to Stroma, the mechanism of "biological" action of its laser would be based on performing sessions lasting 20 seconds, which would cause a slow digestion of iris pigment by macrophages. Thus, although at first the iris would darken, it would become light blue after a month.

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3 COMPLICATIONS OF NEWEYES IN THE SHORT, MEDIUM AND LONG TERM

The technical complications arising from the use of NEWEYES laser are minor, self-limited and resolve on their own within a few hours or days. We refer to fleeting tension elevations stromal micro hemorrhages, mild anterior uveitis, micro iridotomies, transitional anisocorias and short photophobia. All are conveniently prevented and treated by the standard protocol by EYECOS.

Now referring to the pigment deposit or bottom clumps ("6h goop"), which can be observed the first days after the depigmentation sessions, it varies in size depending on the amount of extracted pigment, but inevitably is subjected to a progressive contraction and lessening thus making it appear with a very small sizes after several weeks, or even their complete disappearance.

In the medium term, 1 or 2 years of evolution, the complications observed in patients treated with the NEWEYES laser are simply NULL. No complaints have been reported, such as photophobia, visual disturbances, loss of visual acuity, functional anisocorias, iritis, inflammation or redness nor elevations in intraocular pressure.

The theme of the ocular pressure is really paradoxical. According to the ophthalmologic common sense a slight or moderate elevation of the pressure over time would be expected. In fact many pessimistic colleagues have predicted it. A study by the University of Otago in New Zealand has established the relationship between pigmented eyes and a greater predisposition to glaucoma, but reality reveals just the opposite. Controlled pressures in successive visits, without any topical or oral antihypertensive therapy scheduled, show pressure values that are never higher than the original pressures previous to treatment, and even, on many occasions, somewhat lower.

We have no doubt that the pigment does not block the trabecular meshwork of our young patients (mean age = 33 years). The scattered remnants pass through the trabecular meshwork and ends in the Schlemm's canal after, then being evacuated to the venous bloodstream. On the other hand, much of melanin undergoes agglomeration, forming the "deposit or lower lump" ("6h goop"), which, thanks to the force of gravity stays at that point 6 hours. Over the weeks, this cluster disap-

pears, leaving a small residual mark. Thus, in the worst case that there exists a trabeculum blocking, the latter would only affect a small part of it, an arc of point, the equivalent to the 6 hours staying, the remaining 360° being completely free of pigment. The facts support this hypothesis.

On the other hand, we do not discard a disruptor cycle effect of the NEWEYES laser, which hypothetically could slightly reduce the flow of aqueous humor production, thus explaining the slight decreases in intraocular pressure detected after application of the NEWEYES laser.

The long term complications attributed to NEWEYES relate to the progressive deterioration of the intraocular underlying structures, such as the retina, and more particularly the macula and fovea.

- A clinical review of Eye & Ear Infirmary Massachusetts, has found a higher rates of deaths in systemic metastases in individuals with low pigmented eyes with choroidal melanoma.
- At the University of Detroit a higher prevalence of AMD cases where low iris pigmentation has been found.

If we review the histology of the iris we will see that it is divided into three main layers: the posterior epithelium very pigmented dark brown melanin, the stroma which contains the muscles, vessels and nerves, and the anterior pigmented layer, which always has a tone lighter than the posterior one.

The NEWEYES laser only depigments the outermost surface of the iris, without touching or damaging the posterior stroma or the epithelium. This means that the eye will preserve the main protective barrier against sunlight for the rest of his life, the iris posterior pigment epithelium, present in all human eyes, even in blue eyes by birth.

This fact leads us to believe that the pessimistic predictions that predict a higher rate of AMD (age-related macular degeneration), and a higher rate of metastasis of choroidal melanomas will not be met.

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4 AESTHETIC EXPECTATIONS

The anonymous satisfaction survey conducted recently by EYECOS to all patients treated with NEWEYES has produced very gratifying results in terms of the professionalism of the clinic and the doctor in charge of assistance. The absence of discomfort, problems during the sessions and subsequent complications has been greatly appreciated as well.

By contrast, 50% of the patients manifested not to agree at 100% with the aesthetic result finally obtained. This is because with NEWEYES, THE RESULTING COLOR CANNOT BE CHOSEN. The main advantage of this laser is its naturalness and final discretion. The greenish or blue hue appears progressively after the sessions corresponding to the stromal fibers, and obviously belongs intimately to the anatomy of the subject. The tone will vary by many factors, from dark to extremely bright, but in no case it can be chosen in a Pantone color palette.

Of course, this feature is advised in the Informed Consent that is given to each patient, but those persons concerned, for many years and years, eager for a specific and very determined change in their eyes, are often deaf to the clear warnings from us.

Fortunately, in recent months, we have succeeded in brighten many more the iris treated by using of sodium heparin in the form of eye drops. This substance has a solvent effect in the pigment released into the anterior chamber, and facilitates its evacuation through the trabeculum. It also makes a co adjuvant task for the macrophages, as to the final cleaning of the anterior surface of the iris stroma.

With the use of heparin sodium we have succeeded in changing the natural tendency of patients treated from green towards blue, more or less clear.

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5 FUTURE OF SUPERSPECIALTY: INTRAOCULAR AESTHETICS

Until now, “extraocular” was understood as ocular aesthetics, i.e., the blepharoplasties, largely monopolized by cosmetic surgeons.

The NEWEYES laser paves the way for a new super specialty: INTRAOCULAR AESTHETICS.

Without doubt, the massive and long-awaited possibility to change the eye color foresees a success for this new laser assisted medical technique.

Painstaking and lengthy clinical studies to confirm our working hypothesis as to the safety neweyes will be required. The passing of time and the competition of new ophthalmic professionals will facilitate the improvement of the technique and an increasing degree of confidence in it.

But we do not only anticipate as interesting just the uniform change in eye color. The NEWEYES technology also bring new clinical indications, such as iris tattoo (“eyeTATTOO”), the artistic design of the appearance of the eyes (“neweyes ART”) and iris marking for documentary and identification purposes (“eyeID”).

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6 RESEARCH & DEVELOPMENT

The two main R & D lines of Clinical EYECOS undergo two ambitious projects:

- ➔ The construction of the new version of the NEWEYES laser, most powerful, automated, and with greater versatility and functionality.
- ➔ The development of the new technique “oneday NEWEYES” which enables depigmentation of brown iris with maximum pigmentation without complications and with optimal aesthetic results, and above all, allows to shorten the time required to execute the technique, passing from several sessions in two different weeks to just one or two days in total.

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CONCLUSIONS

- The neweyes laser is a safe, effective and reproducible technique for iris depigmentation.
- The level of early complications is low, self limited and without clinical repercussions.
- The level of complications in the medium term is null.
- The level of long-term complications remains to be seen, but our working hypothesis is that it will be very small, and without sequels.
- The effectiveness of neweyes is 100 %, always reaching depigmentation of the iris.
- The aesthetic result is natural and discrete, varies from green to blue, but it cannot be previously chosen.

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